

# SAPPHIRE PUMP HOMECARE BACKPACK



**Image 1**

- 1. IV Bag area
- 2. IV Bag buckle
- 3. IV Bag support wings
- 4. Sapphire pump pocket
- 5. Administration Set securing shackles/clasps
- 6. Infusion accessories pocket



**Image 2**

- 7. Personal items front pocket
- 8. Transparent Sapphire Window

**Image 3**

- 9. Carrying Handle
- 10. Shoulder straps
- 11. Personal items side pockets

**Do not place personal items in the main pocket of the backpack for risk of interfering with the proper operation of the infusion. Personal items should be kept only in the front pocket [7] or side pockets [11].**

## SET UP

1. Place the IV bag in the IV bag area [1]
2. To secure the IV bag, thread the IV bag buckle [2] through the IV bag handle and clip
3. Fasten the IV bag with the support wings [3], pull the straps to firmly secure and clip
4. Place the Sapphire pump in the designated pocket [4] with the screen facing out of the bag (screen can be seen through the transparent window [8], secure with the strap and clip
5. Secure the administration set with the securing shackles/clips [5]. Select the suitable shackles according to the preferred outlet of the administration set and thread the administration set through those shackles.
6. Close both zippers of the bag, leaving space for the administration set to emerge between them. Do not apply pressure on the set with the zippers
7. The backpack can be carried by hand with the carrying handle [9] or on the back, using the shoulder straps [10]. It can also be placed in an upright position against the wall.

